

Product Datasheet

Human Hp antibody (IgM) ELISA Kit orb407461

Description	Human anti-Helicobacter pylori (Hp) antibody (IgM) ELISA kit
Reactivity	Human
Tested Applications	ELISA
Target	Hp antibody (IgM)
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm