

Product Datasheet

Human APOA2 ELISA Kit orb407481

Description	Human Apolipoprotein A2, apo-A2 ELISA Kit
Reactivity	Human
Range	1 ng /mL-800 ng/mL
Sensitivity	1 ng/mL
Tested Applications	ELISA
Target	APOA2
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P02652
Expiration Date	Please enquire.
Application Notes	• Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm