

Product Datasheet

Human APOC3 ELISA Kit orb407485

Description	Human Apolipoprotein C3,apo-C3ELISA Kit
Reactivity	Human
Range	15.6 ng/ml-1000 ng/ml
Sensitivity	3.9 ng/ml.
Tested Applications	ELISA
Target	APOC3
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm