

Product Datasheet

Human APOM ELISA Kit orb407492

Description	Human Apolipoprotein M(APOM) ELISA kit
Reactivity	Human
Range	1.56 ng/ml-100 ng/ml
Sensitivity	0.39ng/ml
Tested Applications	ELISA
Target	APOM
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm