

Product Datasheet

Human TNFSF13 ELISA Kit orb407503

Description	Human a proliferation inducing ligand(APRIL)ELISA Kit
Reactivity	Human
Range	0.78 ng/ml-50 ng/ml
Sensitivity	0.195 ng/ml
Tested Applications	ELISA
Target	TNFSF13
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, urine• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm