

Product Datasheet

Human HNP1-3 ELISA Kit orb407537

Description	Human neutrophil peptide 1-3,HNP1-3 ELISA Kit
Reactivity	Human
Range	1.56 ng/ml-100 ng/ml
Sensitivity	0.39 ng/ml
Tested Applications	ELISA
Target	HNP1-3
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm