

Product Datasheet

Human TNFAIP6 ELISA Kit orb407557

Description	Human Tumor necrosis factor-inducible gene 6 protein(TNFAIP6) ELISA kit
Reactivity	Human
Tested Applications	ELISA
Target	TNFAIP6
Note	For research use only.
Uniprot ID	P98066
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, cell culture supernates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm