

Product Datasheet

Human TAGLN2 ELISA Kit orb407571

Description	Human Transgelin-2(TAGLN2) ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Sensitivity	0.078 ng/mL
Tested Applications	ELISA
Target	TAGLN2
Note	For research use only.
Sample Types	serum, plasma, cell lysates
Uniprot ID	P37802
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm