

Product Datasheet

Human L-Carnitine ELISA Kit orb407596

| | |
|----------------------------|---|
| Description | Human L-Carnitine ELISA Kit |
| Reactivity | Human |
| Range | 0.312 ng/mL-20 ng/mL |
| Sensitivity | 0.078 ng/mL |
| Tested Applications | ELISA |
| Target | L-Carnitine |
| Note | For research use only. |
| Sample Types | serum, plasma, urine, tissue homogenates |
| Expiration Date | Please enquire. |
| Application Notes | • Sample Type: serum, plasma, urine, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm |