

Product Datasheet

Human 6HMS ELISA Kit orb407598

Description	Human 6-hydroxymelatonin sulfate (6HMS) ELISA kit
Reactivity	Human
Range	10 pg/ml- 400 pg/ml
Sensitivity	10 pg/ml
Tested Applications	ELISA
Target	6HMS
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-3 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450nm