

## Product Datasheet

### Monkey IFN $\alpha$ ELISA Kit orb407648

<b>Description</b>	Monkey Interferon $\alpha$ , IFN- $\alpha$ ELISA Kit
<b>Reactivity</b>	Monkey
<b>Range</b>	15.6 pg/mL-1000 pg/mL
<b>Sensitivity</b>	3.9 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	IFNA1
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, urine, cell culture supernates, tissue homogenates
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	<ul style="list-style-type: none"><li>• Sample Type: serum, urine, cell culture supernates, tissue homogenates</li><li>• Assay Time: 1-5 hours</li><li>• Sample Volume: 50-100 <math>\mu</math>l</li><li>• Detection Wavelength: 450 nm</li></ul>