

## Product Datasheet

### Monkey COR ELISA Kit orb407670

<b>Description</b>	Monkey cortisol(COR) ELISA kit
<b>Reactivity</b>	Monkey
<b>Range</b>	0.049 ng/mL-200 ng/mL
<b>Sensitivity</b>	0.049 ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	COR
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450nm