

Product Datasheet

Mouse PDGFR beta ELISA Kit orb407804

Description	Mouse Beta-type platelet-derived growth factor receptor(PDGFRB) ELISA kit
Reactivity	Mouse
Range	23.44 pg/ml-1500 pg/ml
Sensitivity	5.86 pg/ml.
Tested Applications	ELISA
Target	PDGFR beta
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm