

Product Datasheet

Mouse CA3 ELISA Kit orb407836

Description	Mouse Carbonic anhydrase 3(CA3) ELISA kit
Reactivity	Mouse
Range	4.69 pg/mL-300 pg/mL
Sensitivity	1.17 ng/mL
Tested Applications	ELISA
Target	CA3
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P16015
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm