

Product Datasheet

Mouse Elastin ELISA Kit orb408034

Description	Mouse Elastin(ELN) ELISA kit
Reactivity	Mouse
Range	3.12 ng/ml-200 ng/ml
Sensitivity	0.78 ng/ml
Tested Applications	ELISA
Target	Elastin
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates.• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm