

Product Datasheet

Mouse EGFL6 ELISA Kit orb408056

Description	Mouse Epidermal growth factor-like protein 6(EGFL6) ELISA kit
Reactivity	Mouse
Range	93.75 pg/ml-6000 pg/ml
Sensitivity	23.43 pg/ml
Tested Applications	ELISA
Target	EGFL6
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm