

Product Datasheet

Mouse LGALS1 ELISA Kit orb408116

Description	Mouse Galectin-1(LGALS1) ELISA kit
Reactivity	Mouse
Range	0.312 ng/mL-20 ng/mL
Sensitivity	0.078 ng/mL
Tested Applications	ELISA
Target	LGALS1
Note	For research use only.
Sample Types	serum, plasma,tissue homogenates
Uniprot ID	P16045
Expiration Date	Please enquire.
Application Notes	• Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm