

Product Datasheet

Mouse HAO1 ELISA Kit orb408204

Description	Mouse Hydroxyacid oxidase 1(HAO1) ELISA kit
Reactivity	Mouse
Range	6.25 pg/mL-400 pg/mL
Sensitivity	1.56 pg/mL
Tested Applications	ELISA
Target	HAO1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9WU19
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm