

## Product Datasheet

### Mouse ST2 ELISA Kit orb408240

<b>Description</b>	Mouse Interleukin-1 receptor-like 1(IL1RL1) ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	6.25 pg/mL-400 pg/mL
<b>Sensitivity</b>	1.56 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	IL1RL1
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell culture supernates
<b>Uniprot ID</b>	<b>P14719</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, tissue homogenates, cell culture supernates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm