

Product Datasheet

Mouse WNT5A ELISA Kit orb408544

Description	Mouse Protein Wnt-5a(WNT5A) ELISA kit
Reactivity	Mouse
Range	3.12 pg/ml-200 pg/ml
Sensitivity	0.78 pg/ml
Tested Applications	ELISA
Target	WNT5A
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm