

Product Datasheet

Mouse RSPO2 ELISA Kit orb408590

Description	Mouse R-spondin-2(RSPO2) ELISA kit
Reactivity	Mouse
Range	31.25 pg/mL-2000 pg/mL
Sensitivity	7.8 pg/mL
Tested Applications	ELISA
Target	RSPO2
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8BFU0
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm