

Product Datasheet

Mouse SPARCL1 ELISA Kit orb408652

Description	Mouse SPARC-like protein 1(SPARCL1) ELISA kit
Reactivity	Mouse
Range	78 pg/ml-5000 pg/ml
Sensitivity	19.5 pg/ml
Tested Applications	ELISA
Target	SPARCL1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm