

Product Datasheet

Mouse SDC4 ELISA Kit orb408668

Description	Mouse syndecan 4 (SDC4) ELISA kit
Reactivity	Mouse
Range	3.12 pg/mL-200 pg/mL
Sensitivity	0.78 pg/mL
Tested Applications	ELISA
Target	SDC4
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O35988
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm