

## Product Datasheet

### Mouse Transthyretin ELISA Kit orb408711

<b>Description</b>	Mouse Transthyretin(TTR) ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	3.12 ng/mL-200 ng/mL
<b>Sensitivity</b>	0.78 ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	TTR
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P07309</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma and tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm