

Product Datasheet

Mouse 25(OH)D3/25 HVD3 ELISA Kit orb408776

Description	Mouse 25-hydroxy vitamin D3(25 HVD3) ELISA Kit
Reactivity	Mouse
Range	20 µg/L-100 µg/L
Sensitivity	5µg/L
Tested Applications	ELISA
Target	25(OH)D3/25 HVD3
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm