

Product Datasheet

Mouse HNE ELISA Kit orb408777

Description	Mouse 4-Hydroxynonenal (HNE)ELISA Kit
Reactivity	Mouse
Range	15.6 pmol/mL-1000 pmol/mL
Sensitivity	3.9 pmol/mL
Tested Applications	ELISA
Target	HNE
Note	For research use only.
Sample Types	serum, plasma
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma • Assay Time: 1-5 hours • Sample Volume: 50-100 μ l • Detection Wavelength: 450 nm