

Product Datasheet

Mouse Col III ELISA Kit orb408790

Description	Mouse Collagen Type III, Col III ELISA KIT
Reactivity	Mouse
Range	3.12 pg/mL-200 pg/mL
Sensitivity	0.78 pg/mL
Tested Applications	ELISA
Target	Col III
Note	For research use only.
Sample Types	serum, plasma
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm