

Product Datasheet

Mouse Interferon gamma ELISA Kit orb408817

| | |
|----------------------------|---|
| Description | Mouse C-X-C motif chemokine 11(CXCL11) ELISA kit |
| Reactivity | Mouse |
| Range | 31.25 pg/ml-2000 pg/ml |
| Sensitivity | 7.81 pg/ml |
| Tested Applications | ELISA |
| Target | Interferon gamma |
| Note | For research use only. |
| Expiration Date | Please enquire. |
| Application Notes | <ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, tissue homogenates.• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm |