

Product Datasheet

Mouse KEAP1 ELISA Kit orb408837

Description	Mouse Kelch-like ECH-associated protein 1(KEAP1) ELISA kit
Reactivity	Mouse
Range	4.7 pg/mL-300 pg/mL
Sensitivity	1.17 pg/mL
Tested Applications	ELISA
Target	KEAP1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9Z2X8
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm