

Product Datasheet

Mouse S100B ELISA Kit orb408862

Description	Mouse Protein S100-B(S100B) ELISA kit
Reactivity	Mouse
Range	3.12 pg/ml-200 pg/ml
Sensitivity	0.78 pg/ml
Tested Applications	ELISA
Target	S100B
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm