

Product Datasheet

Mouse IL16 ELISA Kit orb408922

Description	Mouse Interleukin 16,IL-16 ELISA KIT
Reactivity	Mouse
Range	7.8 pg/mL-500 pg/mL
Sensitivity	1.95 pg/mL
Tested Applications	ELISA
Target	IL16
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O54824
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma and tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm