

Product Datasheet

Mouse IL19 ELISA Kit orb408927

Description	Mouse Interleukin-19(IL19) ELISA kit
Reactivity	Mouse
Range	12.5 pg/mL-800 pg/mL
Sensitivity	3.12 pg/mL
Tested Applications	ELISA
Target	IL19
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8CJ70
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm