

Product Datasheet

Mouse IL2 ELISA Kit orb408934

Description	Mouse Interleukin 2,IL-2 ELISA kit
Reactivity	Mouse
Range	15.6 pg/mL-1000 pg/mL
Sensitivity	3.9 pg/mL
Tested Applications	ELISA
Target	IL2
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell culture supernates
Uniprot ID	P04351
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates and cell culture supernates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm