

## Product Datasheet

### Mouse IL27 ELISA Kit orb408940

<b>Description</b>	Mouse Interleukin 27,IL-27 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	31.25 pg/mL-2000 pg/mL
<b>Sensitivity</b>	7.81 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	IL27
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Uniprot ID</b>	<b>Q8K316</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm