

Product Datasheet

Mouse ALT/GPT ELISA Kit orb408971

Description	Mouse Alanine Aminotransferase(ALT/GPT)ELISA Kit
Reactivity	Mouse
Range	0.156 mU/ml-10 mU/ml
Sensitivity	0.039 mU/ml
Tested Applications	ELISA
Target	ALT/GPT
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm