

Product Datasheet

Mouse Myoglobin ELISA Kit orb409122

Description	Mouse Myoglobin,MYO/MB ELISA Kit
Reactivity	Mouse
Range	15.6 ng/mL-1000 ng/mL
Sensitivity	3.9 ng/mL
Tested Applications	ELISA
Target	MB
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P04247
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm