

Product Datasheet

Mouse MMP13 ELISA Kit orb409131

Description	Mouse matrix metalloproteinase 13, MMP-13 ELISA kit
Reactivity	Mouse
Range	93.75 pg/mL-6000 pg/mL
Sensitivity	23.438 pg/mL
Tested Applications	ELISA
Target	MMP13
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P33435
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm