

Product Datasheet

Mouse TIMP2 ELISA Kit orb409138

Description	Mouse tissue inhibitors of metalloproteinase 2, TIMP-2 ELISA kit
Reactivity	Mouse
Range	125 pg/mL-8000 pg/mL
Sensitivity	31.25 pg/mL
Tested Applications	ELISA
Target	TIMP2
Note	For research use only.
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P25785
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, cell culture supernates and tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm