

Product Datasheet

Mouse FGF7 ELISA Kit orb409162

Description	Mouse Keratinocyte Growth Factor(KGF/FGF-7) ELISA Kit
Reactivity	Mouse
Range	1.56 pg/ml-100 pg/ml
Sensitivity	0.39 pg/ml.
Tested Applications	ELISA
Target	FGF7
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm