

## Product Datasheet

### Mouse TNFSF8 ELISA Kit orb409205

<b>Description</b>	Mouse Soluble Cluster of differentiation 30 ligand,sCD30L ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	Request Information
<b>Sensitivity</b>	Request Information
<b>Tested Applications</b>	ELISA
<b>Target</b>	TNFSF8
<b>Note</b>	For research use only.
<b>Uniprot ID</b>	<b>P32972</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm