

Product Datasheet

Mouse TNFSF10 ELISA Kit orb409221

Description	Mouse soluble tumor necrosis factor-related apoptosis inducing ligand,sTRAIL ELISA KIT
Reactivity	Mouse
Range	3.12 pg/ml-200 pg/ml
Sensitivity	0.78 pg/ml
Tested Applications	ELISA
Target	TNFSF10
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm