

Product Datasheet

Mouse Endothelin 1 ELISA Kit orb409270

Description	Mouse Endothelin 1,ET-1 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml-1000 pg/ml
Sensitivity	3.9pg/ml
Tested Applications	ELISA
Target	Endothelin 1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm