

Product Datasheet

Mouse HSP70 ELISA Kit orb409322

Description	Mouse Heat Shock Protein 70,Hsp-70 ELISA Kit
Reactivity	Mouse
Range	0.156 ng/ml-10 ng/ml
Sensitivity	0.039 ng/ml.
Tested Applications	ELISA
Target	HSP70
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm