

Product Datasheet

Mouse proNGF ELISA Kit orb409340

Description	Mouse Pro-Nerve growth factor (proNGF) ELISA kit
Reactivity	Mouse
Range	23.5 pg/mL-1500 pg/mL
Sensitivity	5.8 pg/mL
Tested Applications	ELISA
Target	proNGF
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-3 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450nm