

Product Datasheet

Mouse Eotaxin 3 ELISA Kit orb409365

Description	Mouse Eotaxin 3 ELISA Kit
Reactivity	Mouse
Range	6.25 pg/ml-400 pg/ml
Sensitivity	1.56 pg/ml
Tested Applications	ELISA
Target	Eotaxin 3
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm