

Product Datasheet

Mouse ECP ELISA Kit orb409367

Description	Mouse eosinophil cationic protein,ECP ELISA Kit
Reactivity	Mouse
Range	0.156 ng/mL-10 ng/mL
Sensitivity	0.039 ng/mL
Tested Applications	ELISA
Target	ECP
Note	For research use only.
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm