

Product Datasheet

Mouse AMY2A ELISA Kit orb409490

Description	Mouse pancreatic alpha amylase (AMY2A)ELISA kit
Reactivity	Mouse
Range	9.38 mU/ml-600 mU/ml
Sensitivity	2.3 mU/ml
Tested Applications	ELISA
Target	AMY2A
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm