

Product Datasheet

SQSTM1 (Phospho-Ser28) antibody orb319622

| | |
|-----------------------------|--|
| Description | Rabbit polyclonal antibody to SQSTM1 |
| Species/Host | Rabbit |
| Reactivity | Human |
| Predicted Reactivity | Primate |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Synthetic phospho-peptide corresponding to amino acid residues surrounding Ser28 conjugated to KLH |
| Target | SQSTM1 |
| Form/Appearance | 10 mM HEPES (pH 7.5), 150 mM NaCl, 100 |
| Storage | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| Note | For research use only. |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | Affinity Purified from Pooled Serum |
| MW | 48 kDa |
| Uniprot ID | Q13501 |
| Entrez | 8878 |
| Expiration Date | 12 months from date of receipt. |