

Product Datasheet

Rat CXCR4 ELISA Kit orb410625

Description	Rat CXC-chemokine receptor 4(CXCR4)ELISA Kit
Reactivity	Rat
Range	0.156 ng/ml-10 ng/ml
Sensitivity	0.039 ng/ml
Tested Applications	ELISA
Target	CXCR4
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm