

Product Datasheet

Rat SFTPC ELISA Kit orb410835

Description	Rat Pulmonary Surfactant-associated protein C (SP-C) ELISA kit
Reactivity	Rat
Range	4.69 pg/mL-300 pg/mL
Sensitivity	1.17 pg/mL
Tested Applications	ELISA
Target	SFTPC
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P11685
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm